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ADCC FcγRIIIa(158F) Jurkat Effector Cell Line

ADCC FcγRIIIa(158F) Jurkat 效应细胞

For research use only!

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Reagent	Specification	Manufacturer/Catalogue No.
RPMI 1640	500 mL	Gibco/C11875500BT
Fetal Bovine Serum	500 mL	Cegrogen biotech/A0500-3010
Pen/Strep	100 mL	Thermo/15140-122
Blasticidin	10 mg	Genomeditech/GM-040404-1
Puromycin	25 mg	Genomeditech/GM-040401-1
96 Well Clear V-Bottom Tissue Culture	96-well	Corning/3894
96 well round well culture plate	96-well	NEST/701001
96 well White Flat Bottom Polystyrene Not Treated Microplate	96-well	Corning/3912
H_CLDN18.2(isoform2) CHO-K1 Cell Line		Genomeditech/GM-C05273
Anti-CLDN18.2 hIgG1 Antibody(Zolbetuximab)	/	Genomeditech/GM-34137AB
GMOne-Step Luciferase Reporter Gene Assay Kit	1000T	Genomeditech/GM-040503C

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Equipment	Manufacturer/Catalogue No.
细胞计数仪	ThermoFisher Scientific/Countess 3
酶标仪	Moleculardevices/SpectraMax L

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FcγR	
Cynomolgus_FcRn MDCK Cell Line	H_FCGR1A(CD64) CHO-K1 Cell Line
H_FCGR1A(CD64) HEK-293 Cell Line	H_FCGR2A(CD32A) CHO-K1 Cell Line
H_FCGR2B(CD32B) CHO-K1 Cell Line	H_FCGR3A(CD16a) 158F CHO-K1 Cell Line
H_FCGR3A(CD16a) 158V CHO-K1 Cell Line	H_FCGR3B(CD16b) CHO-K1 Cell Line
H_FcRn CHO-K1 Cell Line	H_FcRn MDCK Cell Line
Mouse_FcRn MDCK Cell Line	
Anti-FcRn hIgG4 Reference Antibody(Rozabio)	Anti-H_FcRn IgG4 Antibody(Rozanolixizumab)
Anti-Mouse CD1632 mIgG2b Antibody(2.4G2)	
ADCCP	
ADCC FcγRIIIa(158V) DDX35TM Jurkat Effector Cell Line	ADCC FcγRIIIa(158V) Jurkat Effector Cell Line
ADCC M_FcγRIV Jurkat Effector Cell Line	ADCP FcγRIIa DDX35TM Jurkat Effector Cell Line
ADCP FcγRIIa Jurkat Effector Cell Line	ADCP FcγRIIa R131 Jurkat Effector Cell Line
ADCP FcγRIIb Jurkat Effector Cell Line	

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