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H_PD-1 Reporter Jurkat Cell Line

H_PD-1 Reporter Jurkat 细胞系

For research use only!

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Reagent	Specification	Manufacturer/Catalogue No.
Puromycin	25 mg	Genomeditech/GM-040401-1
Blasticidin	10 mg	Genomeditech /GM-040404-1
Pen/Strep	100 mL	Thermo/15140-122
Fetal Bovine Serum	500 mL	ExCell/FSP500
RPMI 1640	500 mL	Gibco/C11875500BT
96 Well Clear V-Bottom Tissue Culture	96-well	Corning/3894
96 Well White Polystyrene Microplate	96-well	Corning/3903
Cell Culture Dish	10 cm	NEST/704001
aAPC(OKT3) PDL1 CHO-K1 Cell Line	5E6 cells/mL	Genomeditech/GM-C05269
H_PD-L1 Raji Cell Line	5E6 cells/mL	Genomeditech/GM-C03541
H_CD32B aAPC CHO-K1 Cell Line	5E6 cells/mL	Genomeditech/GM-C25754
OKT3(CD3 ScFv)CHO-K1 Cell Line	5E6 cells/mL	Genomeditech/GM-C19403
金黄色葡萄球菌肠毒素 (SEE 肠毒素)	1 µg	Genomeditech/GM-H23036
Anti-H_CD274(PDL1) hIgG1 Antibody(Atezolizumab)	/	Genomeditech/GM-31740AB
Anti-PD1 hIgG4 Antibody(Pembrolizumab)	/	Genomeditech/GM-52674AB
Anti-PD1 hIgG1 Antibody(Rosnilimab)	/	Genomeditech/GM-87952AB
Anti-PD1 hIgG4 Reference Antibody (Pembio)	/	Genomeditech/GM-87802MAB
GMOne-Step 2.0 Luciferase Reporter Gene Assay Kit	1000T	Genomeditech/GM-040513C

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Equipment	Manufacturer/Catalogue No.
细胞计数仪	ThermoFisher Scientific/Countess 3
酶标仪	Moleculardevices/SpectraMax L

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
Fig 2. Response to Anti-H_CD274(PDL1) hIgG1 Antibody(Atezolizumab). Serial dilutions of the Anti-H_CD274(PDL1) hIgG1 Antibody (Atezolizumab) (Cat. GM-31740AB) was incubated with 1.5E4 cells/well of the aAPC(OKT3) PDL1 CHO-K1 Cell Line (Cat. GM-C05269) in a 96-well plate for 1 hour. Subsequently, H_PD-1 Reporter Jurkat Cell Line (Cat. GM-C07928) with a concentration of 1E5 cells/well was added, and the co-culture proceeded for an additional 16 hours in assay buffer (RPMI 1640 + 1% FBS + 1% P.S). The firefly luciferase activity was measured using the Luciferase Reporter Assay Kit (Genomeditech).

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
Fig 3. Response to Anti-PD1 hIgG4 Antibody(Pembrolizumab). aAPC(OKT3) PDL1 CHO-K1 Cell Line (Cat. GM-C05269) was seeded at a density of 1.5E4 cells/well in a 96-well plate and incubated overnight. The next day, serial dilutions of the Anti-PD1 hIgG4 Antibody(Pembrolizumab) (Cat. GM-52674AB) ,and Human IgG1 Isotype Control were incubated with 1E5 cells/well of the H_PD-1 Reporter Jurkat Cell Line (Cat. GM-C07928) in a 96-well plate for 1 hour, and then added to the pre-seeded cells. The mixture was incubated for an additional 6 hours.

The firefly luciferase activity was measured using the Luciferase Reporter Assay Kit (Genomeditech).

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
Fig 4. Response to Rosnilimab and Isotype. Serial dilutions of the Rosnilimab and Isotype was incubated with $1\text{E}4$ cells/well of the H_CD32B aAPC CHO-K1 Cell Line (Cat. GM-C25754) in a 96-well plate for 1 hour. Subsequently, H_PD-1 Reporter Jurkat Cell Line (Cat. GM-C07928) with a concentration of $1\text{E}5$ cells/well was added, and the co-culture proceeded for an additional 6 hours in assay buffer (RPMI 1640 + 1% FBS + 1% P.S). The firefly luciferase activity was measured using the Luciferase Reporter Assay Kit (Genomeditech).


Fig 5. Response to Rosnilimab and Isotype. Serial dilutions of the Rosnilimab and Isotype was incubated with $1\text{E}4$ cells/well of the OKT3(CD3 ScFv) CHO-K1 Cell Line (Cat. GM-C19403) in a 96-well plate for 1 hour. Subsequently, H_PD-1 Reporter Jurkat Cell Line (Cat. GM-C07928) with a concentration of $1\text{E}5$ cells/well was added, and the co-culture proceeded for an additional 6 hours in assay buffer (RPMI 1640 + 1% FBS + 1% P.S). The firefly luciferase activity was measured using the Luciferase Reporter Assay Kit (Genomeditech).

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
Fig 6. Response to Anti-H_CD274(PDL1) hIgG1 Antibody (Atezolizumab). The Atezolizumab (Cat. GM-31740AB) was serially diluted and incubated with 2E4 cells/well of the H_PD-L1 Raji Cell Line in a 96-well plate for 30 minutes. Meanwhile, 1E5 cells/well of the H_PD-1 Reporter Jurkat Cell Line (Cat. GM-C07928) was pre-incubated with 0.5 ng/well of SEE in a 96-well plate for 30 minutes. Afterward, the two treated cell mixtures were combined in equal volumes and co-incubated in assay buffer (RPMI 1640 + 1% FBS + 1% P.S.) for 16 hours. The firefly luciferase activity was measured using the Luciferase Reporter Assay Kit (Genomeditech). The results indicated a maximum fold of approximately [8.3]. Data are presented based on drug mass concentration.


Fig 7. Response to Anti-PD-1 Human IgG4 Antibody (Pembrolizumab). Having been pre-incubated with 1.5 ng/well of SEE for 30 minutes, 2E4 cells/well of the H_PD-L1 Raji Cell Line (Cat. GM-C03541) were subsequently incubated with serially diluted pembrolizumab for 5 minutes. Upon completion, 1E5 cells/well of the H_PD-1 Reporter Jurkat Cell Line (Cat. GM-C07928) were introduced into the assay buffer, and the reaction was allowed to proceed for 16 hours. The firefly luciferase activity was measured using the Luciferase Reporter Assay Kit (Genomeditech). The results indicated a maximum fold of approximately [2.3]. Data are presented based on drug mass concentration.

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
Fig 8. The passage stability of response to Anti-H_CD274 (PDL1) hIgG1 Antibody (Atezolizumab). Serial dilutions of the Atezolizumab (Cat. GM-31740AB) were incubated with 1.5E4 cells/well of the aAPC (OKT3) PDL1 CHO-K1 Cell Line (Cat. GM-C05269) in a 96-well plate for 1 hour in assay buffer (RPMI 1640 + 1% FBS + 1% P.S). Subsequently, the passage 4, 14, and 24 of H_PD-1 Reporter Jurkat Cell Line (Cat. GM-C07928) at a concentration of 1E5 cells/well was added, and the coculture proceeded for an additional 16 hours. The firefly luciferase activity was measured using the Luciferase Reporter Assay Kit (Genomeditech). Data are shown by drug mass concentration.

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
Fig 9. H_PD-1 Reporter Jurkat Cell Line (Cat. GM-C07928) was determined by flow cytometry using Anti-PD1 hIgG4 Reference Antibody (Pembio) (Cat.GM-87802MAB).

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PD-1:PD-L1(B7-H1):PDL2	
Mouse_PDL1 KO MC38 Cell Line	aAPC(OKT3) PDL1 CHO-K1 Cell Line
H_PDCD1LG2(PDL2) aAPC CHO-K1 Cell Line	Mouse PDL1 aAPC CHO-K1 Cell Line
Mouse_PD-1 Reporter Jurkat Cell Line	Canine_PD-1 HEK-293 Cell Line
Cynomolgus_PD1 CHO-K1 Cell Line	H_CD274(PD-L1) CHO-K1 Cell Line
H_CD274(PD-L1) MC38 Cell Line	H_PDCD1(PD-1) CHO-K1 Cell Line
H_PDCD1LG2(PDL2) CHO-K1 Cell Line	H_PD-L1 HEK-293 Cell Line
H_PDL1 LLC1(mouse_PDL1 KO) Cell Line	H_PDL1 LLC1(mouse_PDL1 KO) Cell Line
H_PDL1 MC38(mouse_PDL1 KO) Cell Line	H_PD-L1 Raji Cell Line
M_PDCD1(PD-1) CHO-K1 Cell Line	
Anti-Canine_PD1 mIgG2a Antibody(4F12-E6)	Anti-H_CD274(PDL1) Antibody(Atezolizumab) hIgG1
Anti-H_PDCD1(PD1) hIgG1 Antibody(Budigalimab)	Anti-H_PDCD1LG2 mIgG1 Antibody(3G2)
Anti-mouse PD1 RIgG2a Antibody(RMP1-14)	Anti-mouse PD-L1 mIgG1 Antibody(10F.9G2)
Anti-Mouse_PD1 mIgG1 Antibody(29F.1A12)	Anti-mouse_PD1 mIgG1 Antibody(RMP1-14)
Anti-PD1 hIgG4 Antibody(Pembrolizumab)	Anti-PD1 hIgG4 Reference Antibody (Nivbio)
Anti-PD1 hIgG4 Reference Antibody (Pembio)	Anti-PD1 hIgG4 Reference Antibody (Sintbio)
Anti-PD-1 hIgG4 Reference Antibody (Torbio)	Anti-PD1 hIgG4 Reference Antibody(Cambio)
Anti-PD-1 hIgG4 Reference Antibody(Tislbio)	Anti-PD-L1 hIgG1 Reference Antibody(Avebio)
Anti-PDL1 hIgG4 Reference Antibody(Adebio)	Anti-PD-L2 hIgG1 Antibody(Hz25G4-1.1)
Biotinylated Human PD1 Protein; His-Avi Tag	Biotinylated Human PDL1 Protein; His-Avi Tag
Canine PD1 Protein; hFc Tag	Cynomolgus PDL1 Protein; His Tag
Human PD1 Protein; His Tag	Human PDL1 Protein; His Tag
Mouse PDL1 Protein; His Tag	

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