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ADCC FcγRIIIa(158V) Jurkat Effector Cell Line

ADCC FcγRIIIa(158V) Jurkat 效应细胞

For research use only!

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Reagent	Specification	Manufacturer/Catalogue No.
Puromycin	25mg	Genomeditech/GM-040401-1
Blasticidin	10mg	Genomeditech/GM-040404-1
Pen/Strep	100 mL	Thermo/15140-122
Fetal Bovine Serum	500 mL	ExCell/FSP500
RPMI-1640	500 mL	Viva Cell/C3010-0500
96 Well Clear V-Bottom Tissue Culture	96-well	Corning/3894
96 well round well culture plate	96-well	NEST/701001
96 well White Flat Bottom Polystyrene	96-well	Corning/3912
Not Treated Microplate		
GMOne-Step 2.0 Luciferase Reporter Gene Assay Kit	1000T	Genomeditech/GM-040513C
Anti-CLDN18.2 hIgG1 Antibody	/	Genomeditech/GM-34137AB
Anti-CD20 (Rituximab)	/	Expression
PE anti-human CD16 Antibody	25 tests	biolegend/302007

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Equipment	Manufacturer/Catalogue No.
细胞计数仪	ThermoFisher Scientific/Countess 3
酶标仪	Moleculardevices/SpectraMax L

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
Fig 2. Response to Anti-CLDN18.2 hIgG1 Antibody. Serial dilutions of the Anti-CLDN18.2 hIgG1 Antibody(Zolbetuximab) (Cat. GM-34137AB) and 1.5E5 cells/well of the ADCC FcγRIIIa(158V) Jurkat Effector Cell Line (Cat. GM-C05619) were added to 1.5E4 cells/well of H_CLDN18.2 CHO-K1 cell line (Cat. GM-C05273) for 6 hours. Firefly luciferase activity is then measured using the Luciferase Reporter Assay Kit (Genomeditech). The maximum induction fold was approximately[49.8]. Data are shown by drug mass concentration.

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
Fig 3. Response to Anti-CD20(Rituximab). Serial dilutions of the Anti-CD20(Rituximab) and 1.5E5 cells/well of the ADCC Fc γ RIIIa(158V) Jurkat Effector Cell Line (Cat. GM-C05619) were added to 2.5E4 cells/well of Raji cell line (Cat. GM-C19100) for 6 hours. Firefly luciferase activity is then measured using the Luciferase Reporter Assay Kit (Genomeditech). The maximum induction fold was approximately[45.9]. Data are shown by drug molar concentration.

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
Fig 4. Response to Anti-CLDN18.2 hIgG1 Antibody. Serial dilutions of the Anti-CLDN18.2 hIgG1 Antibody(Zolbetuximab) (Cat. GM-34137AB) and 1.5E5 cells/well of the passage 5, 16 and 28 of the ADCC FcγRIIIa(158V) Jurkat Effector Cell Line (Cat. GM-C05619) were added to 1.5E4 cells/well of H₂CLDN18.2 CHO-K1 cell line (Cat. GM-C05273) for 6 hours. Firefly luciferase activity is then measured using the Luciferase Reporter Assay Kit (Genomeditech). Data are shown by drug molar concentration.

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ADCC FcγRIIIa(158V) Jurkat Effector Cell Line+175D10	0 μg/mL	10 μg/mL	2.38 pg/mL
	288	8716	361
ADCC FcγRIIIa(158V) Jurkat Effector Cell Line+Enhanced 175D10	0 μg/mL	10 μg/mL	2.38 pg/mL
	285	10095	289
ADCC FcγRIIIa(158V) Jurkat Effector Cell Line+Antibody-1	0 μg/mL	10 μg/mL	2.38 pg/mL
	310	7018	205
ADCC FcγRIIIa(158V) Jurkat Effector Cell Line+Antibody-2	0 μg/mL	10 μg/mL	2.38 pg/mL
	299	7531	296
ADCC FcγRIIIa(158V) Jurkat Effector Cell Line+Igg	0 μg/mL	10 μg/mL	2.38 pg/mL
	265	533	375

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
Fig 6. Raji cell line (Cat. GM-C19100) was determined by flow cytometry using Anti-H_MS4A1(CD20) hIgG1 Antibody(Ocrelizumab) (Cat. GM-27200AB).

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
Fig 7. ADCC FcγRIIIa(158V) Jurkat Effector Cell Line (Cat. GM-C05619) was determined by flow cytometry using PE anti-human CD16 Antibody (BioLegend/302007).

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
Fig 8. Anti-VEGF-mediated ADCC assays on CHO-K1 cells. The ADCC FcγRIIIa(158V) Jurkat Effector Cell Line(Cat. GM-C05619) at a concentration of 1E5 cells/well was co-cultured with Membrane Bound H₂VEGF165 CHO-K1 Cell Line (Cat. GM-C40448) at a concentration of 1.5E4 cells/well, in the presence of serial dilutions of the Anti-VEGF hIgG1 Reference Antibody (Bevbio) (Cat.GM-87758MAB) for 6 hours (96-well format). The firefly luciferase activity was measured using the Luciferase Reporter Assay Kit (Genomeditech). Data are shown by drug mass concentration.

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FcγR	
Cynomolgus_FcRn MDCK Cell Line	H_FCGR1A(CD64) CHO-K1 Cell Line
H_FCGR1A(CD64) HEK-293 Cell Line	H_FCGR2A(CD32A) CHO-K1 Cell Line
H_FCGR2B(CD32B) CHO-K1 Cell Line	H_FCGR3A(CD16a) 158F CHO-K1 Cell Line
H_FCGR3A(CD16a) 158V CHO-K1 Cell Line	H_FCGR3B(CD16b) CHO-K1 Cell Line
H_FcRn CHO-K1 Cell Line	H_FcRn MDCK Cell Line
Mouse_FcRn MDCK Cell Line	
Anti-FcRn hIgG4 Reference Antibody(Rozabio)	Anti-H_FcRn IgG4 Antibody(Rozanolixizumab)
Anti-Mouse CD1632 mIgG2b Antibody(2.4G2)	
ADCCP	
ADCC FcγRIIIa(158F) Jurkat Effector Cell Line	ADCC FcγRIIIa(158V) DDX35TM Jurkat Effector Cell Line
ADCC M_FcγRIV Jurkat Effector Cell Line	ADCP FcγRIIa DDX35TM Jurkat Effector Cell Line
ADCP FcγRIIa Jurkat Effector Cell Line	ADCP FcγRIIa R131 Jurkat Effector Cell Line
ADCP FcγRIIb Jurkat Effector Cell Line	

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