Genomeditech (Shanghai) Co.,Ltd. Order: 021-68455258/50432826/50432825

Toll-free: 400 627 9288

Email: service@genomeditech.com

Anti-TREM1 hlgG1 Antibody

Product information

GM-26835AB-10 10 μg GM-26835AB-100 100 μg GM-26835AB-1000 1 mg

Antibody Information

Species Reactivity Human;

Clone

Source/Isotype Monoclonal human IgG1, κ

Application Flow cytometry; Activation assay;

Specificity Detects TREM1

Gene TREM1

Other Names CD354, TREM-1
Gene ID 54210 (human);

Background This gene encodes a receptor belonging to the lg superfamily that is

expressed on myeloid cells. This protein amplifies neutrophil and monocyte-mediated inflammatory responses triggered by bacterial and fungal infections by stimulating release of pro-inflammatory chemokines and cytokines, as well as increased surface expression of cell activation markers. Alternatively spliced transcript variants encoding different

isoforms have been noted for this gene.

Storage Storage Store at 2-8°C short term (1-2 weeks). Store at ≤ -20°C long term. Avoid

repeated freeze-thaw.

Formulation Phosphate-buffered solution, pH 7.2.

Endotoxin < 1 EU/mg, determined by LAL gel clotting assay

Version:3.2 Revision



Genomeditech (Shanghai) Co.,Ltd. Order: 021-68455258/50432826/50432825


Toll-free: 400 627 9288

Email: service@genomeditech.com

Data Examples

Flow cytometry

H_TREM1 CHO-K1 Cell Line (Catalog # GM-C13367) was stained with Anti-TREM1 hlgG1 Antibody (Catalog # GM-26835AB) or isotype control antibody, followed by anti-Human lgG APC-conjugated Secondary Antibody.




SampleID	Geometric Mean : B530-H
CHO-K1 anti-H_TREM1+FITC-2nd Ab	895
CHO-K1 H_TREM1 H_IgG+FITC-2nd Ab	1419
CHO-K1 H_TREM1 anti-H_TREM1+FITC-2nd Ab	384726

Fig. FACS

Flow cytometry

H_TREM1 Reporter Jurkat Cell Line (Catalog # GM-C15720) was stained with Anti-TREM1 hlgG1 Antibody (Catalog # GM-26835AB) or isotype control antibody, followed by anti-Human lgG APC-conjugated Secondary Antibody.



	SampleID	Geometric Mean : BL1-H	
	Jurkat anti-H_TREM1+FITC-2nd Ab	110	
	H_TREM1 Reporter Jurkat H_IgG+FITC-2nd Ab	134	
	H_TREM1 Reporter Jurkat anti-H_TREM1+FITC-2nd Ab	96771	

Fig. FACS



Genomeditech (Shanghai) Co.,Ltd. Order: 021-68455258/50432826/50432825

Toll-free: 400 627 9288

Email: service@genomeditech.com

Activation assay

Anti-TREM1 hlgG1 Antibody (Catalog # GM-26835AB) activates the H_TREM1 Reporter Jurkat Cell Line (Catalog # GM-C15720) after coating, leading to induced luminescence. The EC50 for this effect is $1.808 \ \mu g/mL$.




Fig. assay