

Anti-H_HAVCR2(TIM3) hIgG4 Antibody(Sabatolimab)

Product information

GM-27128AB-10	10 µg
GM-27128AB-100	100 µg
GM-27128AB-1000	1 mg

Antibody Information

Species Reactivity	Human
Clone	Sabatolimab
Source/Isotype	Monoclonal human IgG4, κ
Application	Flow cytometry
Other Names	CD366, HAVcr-2, KIM-3, SPTCL, TIMD-3, TIMD3, Tim-3
Gene ID	NM_032782.5(human)
Storage	Store at 2-8°C short term (1-2 weeks).Store at ≤ -20°C long term. Avoid repeated freeze-thaw.
Formulation	Phosphate-buffered solution, pH 7.2.
Endotoxin	< 1 EU/mg, determined by LAL gel clotting assay

Data Examples

Flow cytometry H_HAVCR2(TIM3) CHO-K1 Cell Line(Catalog # GM-C19183) was stained with Anti-H_HAVCR2(TIM3) hIgG4 Antibody(Sabatolimab) or isotype control antibody.




Fig. FACS

Version:0.3.2.240521