

## Anti-H\_PVRIG hlgG1 Antibody

### Product information


GM-27659AB-10	10 µg
GM-27659AB-100	100 µg
GM-27659AB-1000	1 mg

### Antibody Information

Species Reactivity	Human
Source/Isotype	Monoclonal human IgG1, κ
Application	Flow cytometry;
Other Names	C7orf15, CD112R
Gene ID	Q6DKI7-1(human)
Storage	Store at 2-8°C short term (1-2 weeks).Store at ≤ -20°C long term. Avoid repeated freeze-thaw.
Formulation	Phosphate-buffered solution, pH 7.2.
Endotoxin	< 1 EU/mg, determined by LAL gel clotting assay

### Data Examples

Flow cytometry H\_CD112R(PVRIG) CHO-K1 Cell Line (Catalog # GM-C18950) was stained with Anti-H\_PVRIG hlgG1 Antibody or isotype control antibody.



SampleID	Geometric Mean : RL1-H
CHO-K1 anti-H_CD112R+APC-2nd Ab	116
CHO-K1 H_CD112R H_IgG+APC-2nd Ab	65.9
CHO-K1 H_CD112R anti-H_CD112R+APC-2nd Ab	1853

Fig. FACS

Version:0.3.2.240521