

## Anti-SARS-CoV-2 hIgG1 Antibody (Imdevimab)

### Product information

GM-28192AB-10	10 µg
GM-28192AB-100	100 µg
GM-28192AB-1000	1 mg

### Antibody Information

Species Reactivity	Human;
Clone	Imdevimab
Source/Isotype	Monoclonal human IgG1, λ
Application	Block assay;
Storage	Store at 2-8°C short term (1-2 weeks).Store at ≤ -20°C long term. Avoid repeated freeze-thaw.
Formulation	Phosphate-buffered solution, pH 7.2.
Endotoxin	< 1 EU/mg, determined by LAL gel clotting assay

**Data Examples**

Block assay

Neutralization of wild-type SARS-CoV-2 pseudoviruses




Fig. assay

Block assay

Neutralization of SARS-CoV-2 (Delta) pseudoviruses




Fig. assay

Block assay

Neutralization of SARS-CoV-2 (Delta+) pseudoviruses




Fig. assay