

## Anti-CTGF hIgG1 Antibody(Pamrevlumab)

### Product information


GM-51031AB-10	10 µg
GM-51031AB-100	100 µg
GM-51031AB-1000	1 mg

### Antibody Information

Species Reactivity	Human;
Clone	Pamrevlumab
Source/Isotype	Monoclonal human IgG1, κ
Application	Flow cytometry;
Other Names	CTGF, HCS24, IBP-8, IGFBP8, NOV2
Gene ID	1490(human)
Storage	Store at 2-8°C short term (1-2 weeks).Store at ≤ -20°C long term. Avoid repeated freeze-thaw.
Formulation	Phosphate-buffered solution, pH 7.2.
Endotoxin	< 1 EU/mg, determined by LAL gel clotting assay

### Data Examples

Flow cytometry H\_CTGF CHO-K1 Cell Line was stained with Anti-CTGF hIgG1 Antibody(Pamrevlumab) or isotype control antibody.



SampleID	Geometric Mean : FL11-H
CHO-K1 anti-CTGF+APC-2nd Ab	1395
CHO-K1 H_CTGF H_IgG+APC-2nd Ab	1676
CHO-K1 H_CTGF anti-CTGF+APC-2nd Ab	90751

Fig. FACS

Version:0.3.2.240521