

Anti-DLL4 hIgG2 Antibody(Navicixizumab)

Product information

GM-51815AB-10	10 µg
GM-51815AB-100	100 µg
GM-51815AB-1000	1 mg

Antibody Information

Species Reactivity	Human;
Clone	Navicixizumab
Source/Isotype	Monoclonal human IgG2, κ
Application	Flow cytometry;
Other Names	AOS6, delta4, hdelta2
Gene ID	Q9NR61-1 (human);
Storage	Store at 2-8°C short term (1-2 weeks).Store at ≤ -20°C long term. Avoid repeated freeze-thaw.
Formulation	Phosphate-buffered solution, pH 7.2.
Endotoxin	< 1 EU/mg, determined by LAL gel clotting assay

Data Examples

Flow cytometry H_DLL4 CHO-K1 Cell Line (Catalog # GM-C19190) was stained with Anti-DLL4 hIgG2 Antibody(Navicixizumab) or isotype control antibody.




Fig. FACS

Version:0.3.2.240521