Genomeditech (Shanghai) Co.,Ltd. Order: 021-68455258/50432826/50432825

Toll-free: 400 627 9288

Email: service@genomeditech.com

Anti-HER2-DXD ADC(Dar8)[Trastuzumab Deruxtecan]

Product Information

Product Name Anti-HER2-DXD ADC(Dar8)[Trastuzumab Deruxtecan]

Storage temp. Store at +4°C short term. Store at -80°C for 36 months, avoiding

freeze/thaw cycles.

Catalog# / Size GM-81145AB-100 / 100 μg

GM-81145AB-1mg / 1 mg

Antibody Information

Expression System CHO

Aggregation < 5% as determined by SEC-HPLC

Purity > 95% as determined by SDS-PAGE

Drug-to-Antibody Ratio (DAR) 8.0

Endotoxin <1 EU/mg, determined by LAL gel clotting assay

Sterility 0.22 µm Filtered

Target HER2

Clone Trastuzumab Deruxtecan
Alternative Names DS-8201; DS-8201a

Source/Isotype Monoclonal Human IgG1 (KEEM),Kappa
Application Positive control of Cytotoxicity Assay

Description Trastuzumab deruxtecan (DS-8201) is an antibody-drug conjugate

composed or ADC of an anti-HER2 (human epidermal growth factor receptor 2) antibody, a cleavable tetrapeptide-based linker,

and a cytotoxic topoisomerase I inhibitor.


Formulation 20 mM Histidine, pH 6.0

Genomeditech (Shanghai) Co.,Ltd. Order: 021-68455258/50432826/50432825

Toll-free: 400 627 9288


Email: service@genomeditech.com

SDS PAGE




After SDS-PAGE under reducing conditions, the gel was stained with Coomassie Blue overnight. The purity of the protein is greater than 95%.

SEC-HPLC



The purity of this product is greater than 95% verified by SEC-HPLC.

RP-HPLC



The drug-to-antibody ratio (DAR) of this product is 7.44