Genomeditech (Shanghai) Co.,Ltd. Order: 021-68455258/50432826/50432825

Toll-free: 400 627 9288

Email: service@genomeditech.com

## Anti-ITGB6 hlgG1 Antibody(SGN-B6A)

## **Product information**

GM-81622AB-10 10 μg GM-81622AB-100 100 μg GM-81622AB-1000 1 mg

## **Antibody Information**

Species Reactivity Human
Clone SGN-B6A

Source/Isotype Monoclonal human IgG1/κ

Application Flow cytometry
Specificity Detects ITGB6

Gene ITGB6

Other Names Integrin beta-6
Gene ID 3694(Human)

Background ITGB6 is the  $\beta$ 6 subunit of the integrin  $\alpha v \beta$ 6, which is a crucial subunit

controlling the expression of the entire heterodimer  $\alpha\nu\beta6$ . (ITGB6) is the limiting subunit for the formation of the  $\alpha\nu\beta6$  integrin, as it can only bind with  $\alpha\nu$ . Therefore, the expression of the  $\alpha\nu\beta6$  integrin depends on the transcriptional activity of the  $\beta6$  integrin gene, which is only expressed in epithelial cells. The integrin  $\alpha\nu\beta6$  is exclusively expressed on epithelial cells and is expressed at low levels in most normal healthy tissues. However, it is frequently upregulated in diseases such as cancer and fibrosis, with an estimated one-third of all solid tumors showing upregulation, which is associated with disease progression and poor prognosis. Multiple studies have shown that the severity of many cancers, including lung, breast, colon, gastric, and ovarian cancer, is related to the expression of  $\alpha\nu\beta6$ . In summary,  $\alpha\nu\beta6$  plays a regulatory role in epithelial remodeling during development, tissue repair, and tumor formation.

Storage Store at 2-8°C short term (1-2 weeks). Store at ≤ -20°C long term. Avoid

repeated freeze-thaw.

Formulation Phosphate-buffered solution, pH 7.2.

Endotoxin < 1 EU/mg, determined by LAL gel clotting assay

Version: 3.1 Revision Date: 12/25/2023



Genomeditech (Shanghai) Co.,Ltd.
Order: 021-68455258/50432826/50432825


Toll-free: 400 627 9288

Email: service@genomeditech.com

## **Data Examples**

Flow cytometry

The recommended usage range is 0.5-4  $\mu g$  per test. H\_ITGB6 CHO-K1 Cell Line(Catalog # GM-C29044) was stained with Anti-ITGB6 hlgG1 Antibody(SGN-B6A) (Catalog # GM-81622AB) or isotype control antibody, followed by anti-Human lgG APC-conjugated Secondary Antibody .



SampleID	Geometric Mean : FL11-H
CHO-K1 anti-ITGB6+APC-2nd Ab	1977
CHO-K1 H_ITGB6 H_IgG+APC-2nd Ab	1733
CHO-K1 H_ITGB6 anti-ITGB6+APC-2nd Ab	92617

Fig. FACS