Toll-free: 400 627 9288

Email: service@genomeditech.com

Anti-B7-H6 hlgG1 Antibody(BI-765049)

Product information

GM-86750AB-10 10 μg GM-86750AB-100 100 μg GM-86750AB-1000 1 mg

Antibody Information

Species Reactivity Human;
Clone BI-765049

Source/Isotype Monoclonal human IgG1, κ

Application Flow cytometry; Bioactivity-ELISA

Specificity Detects B7-H6

Gene B7-H6

Other Names NCR3LG1 , B7H6, B7-H6, DKFZp686O24166

Gene ID 374383 (human); 102121659 (Cynomolgus)

Background B7-H6 is a type I transmembrane protein belonging to the B7 family,

characterized by two immunoglobulin domains. It also serves as a ligand for the NK cell-activating receptor NKp30. While B7-H6 is prominently expressed on various tumor cells, no B7-H6 mRNA has been detected in normal human tissues or healthy peripheral blood mononuclear cells. In an inflammatory milieu, certain pro-inflammatory cytokines like IL-1β and TNF can induce the upregulation of B7-H6 expression in CD14+CD16+ monocytes and neutrophils. The tumor expression spectrum of B7-H6 is extensive, encompassing leukemia, lymphoma, colorectal cancer, non-small cell lung cancer, breast cancer, ovarian cancer, stomach cancer, and

liver cancer.

Storage Store at 2-8°C short term (1-2 weeks). Store at ≤ -20°C long term. Avoid

repeated freeze-thaw.

Formulation Phosphate-buffered solution, pH 7.2.

Endotoxin < 1 EU/mg, determined by LAL gel clotting assay

Version:3.4


Toll-free: 400 627 9288

Email: service@genomeditech.com

Data Examples

Flow cytometry

H_NCR3LG1(B7-H6) CHO-K1 Cell Line (Catalog #GM-C36446)was stained with Anti-B7-H6 hlgG1 Antibody(BI-765049) (Catalog # GM-86750AB) or isotype control antibody, followed by anti-Human lgG APC-conjugated Secondary Antibody.




П	SampleID	Geometric Mean : FL11-H
I	HEK-293 anti-B7-H6+APC-2nd Ab	2944
	HEK-293 H_NCR3LG1(B7-H6) H_IgG+APC-2nd Ab	3376
ı	HEK-293 H_NCR3LG1(B7-H6) anti-B7-H6+APC-2nd Ab	1.42E6

Fig. FACS

Flow cytometry

Cynomolgus_NCR3LG1(B7-H6) CHO-K1 Cell Line (Catalog # GM-C36448)was stained with Anti-B7-H6 hlgG1 Antibody(BI-765049) (Catalog # GM-86750AB) or isotype control antibody, followed by anti-Human lgG APC-conjugated Secondary Antibody.



SampleID	Geometric Mean : FL11-H	
CHO-K1 anti-B7-H6+APC-2nd Ab	971	
CHO-K1 Cyno_NCR3LG1(B7-H6) H_IgG+APC-2nd Ab	1007	
CHO-K1 Cyno_NCR3LG1(B7-H6) anti-B7-H6+APC-2nd Ab	289579	

Fig. FACS

Toll-free: 400 627 9288

Email: service@genomeditech.com

Bioactivity-ELISA

Human NCR3LG1(B7-H6) Protein; His Tag (Catalog # GM-87646RP) was immobilized at 1 μ g/ml (100 μ L/well). Increasing concentrations of Anti-B7-H6 hlgG1 Antibody(BI-765049) (Catalog # GM-86750AB) were added.

Bioactivity-ELISA

0.1 µg Human NCR3LG1(B7-H6) Protein; His Tag of per well




Fig. assay

Bioactivity-ELISA

Cynomolgus NCR3LG1(B7-H6) Protein; His Tag (Catalog # GM-87649RP) was immobilized at 1 μ g/ml (100 μ L/well). Increasing concentrations of Anti-B7-H6 hlgG1 Antibody(BI-765049) (Catalog # GM-86750AB) were added.

Bioactivity-ELISA

0.1 μg Cynomolgus NCR3LG1(B7-H6) Protein; His Tag of per well




Fig. assay



Toll-free: 400 627 9288

Email: service@genomeditech.com

Bioactivity-ELISA

Biotinylated Human NCR3LG1(B7-H6) Protein; His-Avi Tag (Catalog # GM-87647RP) was immobilized at 1 μ g/ml (100 μ L/well) on streptavidin precoated. Increasing concentrations of Anti-B7-H6 hlgG1 Antibody(BI-765049) (Catalog # GM-86750AB) were added.

Bioactivity-ELISA

0.1 μg Biotinylated Human NCR3LG1(B7-H6) Protein; His-Avi Tag of per well




Fig. assay

Bioactivity-ELISA

Biotinylated Cynomolgus NCR3LG1(B7-H6) Protein; His-Avi Tag (Catalog # GM-87650RP) was immobilized at 1 μ g/ml (100 μ L/well) on streptavidin precoated. Increasing concentrations of Anti-B7-H6 hlgG1 Antibody(BI-765049) (Catalog # GM-86750AB) were added.

Bioactivity-ELISA

0.1 μg Biotinylated Cynomolgus NCR3LG1(B7-H6) Protein; His-Avi Tag of per well




Fig. assay