Genomeditech (Shanghai) Co.,Ltd. Order: 021-68455258/50432826/50432825

Toll-free: 400 627 9288

Email: service@genomeditech.com

Anti-DLL3 hlgG1 Reference Antibody (Rovabio)

Product Information

Product Name Anti-DLL3 hlgG1 Reference Antibody (Rovabio)

Storage temp. Store at 2-8°C short term (1-2 weeks). Store at ≤ -20°C long term. Avoid

repeated freeze-thaw.

Catalog# / Size GM-87163MAB-1mg / 1 mg

> GM-87163MAB-5mg / 5 mg GM-87163MAB-25mg / 25 mg GM-87163MAB-50mg / 50 mg GM-87163MAB-100mg / 100 mg

Antibody Information

Expression System CHO

Aggregation < 5% as determined by SEC-HPLC **Purity** > 95% as determined by SDS-PAGE

Endotoxin < 1 EU/mg, determined by LAL gel clotting assay

Sterility 0.2 µm Filtered

Target DLL3

Clone Rovalpituzumab

Alternative Names SCD01

Source/Isotype Human IgG1(KDEL, delta 447K), Kappa

Application Flow cytometry

Description Delta-like 3 (Drosophila), also known as DLL3, is a protein which in

> humans is encoded by the DLL3 gene. This gene encodes a member of the delta protein ligand family. This family functions as Notch ligands that are characterized by a DSL domain, EGF repeats, and a transmembrane domain. Expression of DLL3 is highest in fetal brain. It plays a key role in somitogenesis within the Paraxial mesoderm. Mutations in this gene cause the autosomal recessive genetic disorder Jarcho-Levin syndrome. Expression of the gene occurs in Neuroendocrine tumors, which has been targeted as a potential pathway for treatment. An experimental drug, rovalpituzumab tesirine, targets DLL3 as a possible treatment for lung

cancer.

Formulation phosphate-buffered solution, pH 7.4.

Version:3.2


Genomeditech (Shanghai) Co.,Ltd.
Order: 021-68455258/50432826/50432825

Toll-free: 400 627 9288

Fmail: service@genomeditech.com


Data Examples

SDS-PAGE



On SDS-PAGE under reducing (R)/non-reducing(N-R) condition. The gel was stained overnight with Coomassie Blue. The purity of the protein is greater than 95%.

SEC-HPLC



The purity of this product is more than 95% verified by SEC-HPLC

Genomeditech (Shanghai) Co.,Ltd. Order: 021-68455258/50432826/50432825


Toll-free: 400 627 9288

Fmail: service@genomeditech.com

Data Examples

Flow cytometry

H_DLL3 HEK-293 Cell Line (Catalog # GM-C12896) was stained with Anti-DLL3 hlgG1 Reference Antibody (Rovabio) (Catalog # GM-87163MAB) or isotype control antibody, followed by anti-Human lgG APC-conjugated Secondary Antibody.



SampleID	Geometric Mean : FL11-H
HEK-293 anti-DLL3+APC-2nd Ab	1092
HEK-293 H_DLL3 H_lgG+APC-2nd Ab	1096
HEK-293 H_DLL3 anti-DLL3+APC-2nd Ab	53296

Fig. FACS