Genomeditech (Shanghai) Co.,Ltd. Order: 021-68455258/50432826/50432825

Toll-free: 400 627 9288

Email: service@genomeditech.com

Anti-STEAP1 hlgG1 Reference Antibody (Vandbio)

Product Information

Product Name Anti-STEAP1 hlgG1 Reference Antibody (Vandbio)

Storage temp. Store at 2-8°C short term (1-2 weeks). Store at ≤ -20°C long term. Avoid

repeated freeze-thaw.

Catalog# / Size GM-87355MAB-1mg / 1 mg

GM-87355MAB-5mg / 5 mg GM-87355MAB-25mg / 25 mg GM-87355MAB-50mg / 50 mg GM-87355MAB-100mg / 100 mg

Antibody Information

Expression System CHO

Aggregation < 5% as determined by SEC-HPLC

Purity > 95% as determined by SDS-PAGE

Endotoxin < 1 EU/mg, determined by LAL gel clotting assay

Sterility 0.2 µm Filtered

Target STEAP1

Clone Vandortuzumab

Alternative Names PRSS24, STEAP

Source/Isotype Human IgG1 (KEEM), Kappa

Application /

Description The STEAP1 (Six-Transmembrane Epithelial Antigen of Prostate 1) gene

is a gene that is highly expressed in prostate tissue and belongs to a family of six-transmembrane proteins. The main functions of STEAP1 are associated with the transport across the cell membrane, signal

transduction, and interactions between cells.

Formulation phosphate-buffered solution, pH 7.4.


Genomeditech (Shanghai) Co.,Ltd. Order: 021-68455258/50432826/50432825

Toll-free: 400 627 9288

Email: service@genomeditech.com


Data Examples

SDS-PAGE



On SDS-PAGE under reducing (R)/non-reducing(N-R) condition. The gel was stained overnight with Coomassie Blue. The purity of the protein is greater than 95%.

SEC-HPLC



The purity of this product is more than 95% verified by SEC-HPLC