Genomeditech (Shanghai) Co.,Ltd.
Order: 021-68455258/50432826/50432825

Toll-free: 400 627 9288

Email: service@genomeditech.com

# Human PD1 Protein; His Tag

#### **Product Information**

Product Name Human PD1 Protein; His Tag

Storage temp Store at  $\leq$  -70°C, stable for 6 months after receipt.

Recommend to aliquot the protein into smaller quantities for

optimal storage. Please minimize freeze-thaw cycles.

Catalog# / Size GM-87593RP-100 / 100 μg

GM-87593RP-1000 / 1 mg

#### **Protein Information**

Alternative Names PDCD1,PD1,CD279,SLEB2

Source Human PD1 Protein; His Tag (GM-87593RP) is expressed from human 293 cells

(HEK-293). It contains AA Leu 25 - Gln 167 (Accession # Q15116-1).

This protein carries a His tag at the C-terminus.

Purity > 95% as determined by SDS-PAGE

Endotoxin < 1 EU/μg, determined by LAL gel clotting assay

Predicted Mol Mass 16.4 KDa

**Formulation** Supplied as a 0.2 μm filtered solution of PBS, pH7.4.

**Description** PD-1 (Programmed cell death protein 1), also known as PD1 protein, is an

immune regulatory receptor located on the surface of human cells. PD-1 is an immune checkpoint protein that regulates the activation and function of T cells by interacting with its ligands, PD-L1 (Programmed Cell Death Ligand 1) or PD-L2. PD-1 is primarily expressed on activated T cells and B cells, playing a role in

suppressing immune responses and maintaining immune system balance.

Research has shown that PD-1 protein plays an important role in the tumor microenvironment, where tumor cells can evade attacks from the immune system by binding to PD-L1. Therefore, the interaction between PD-1 and PD-L1 has become a critical target in tumor immunotherapy. Inhibiting the interaction between PD-1 and PD-L1 can restore the immune cells' ability to recognize and


attack tumors, providing a new approach for tumor immunotherapy.

Genomeditech (Shanghai) Co.,Ltd. Order: 021-68455258/50432826/50432825

Toll-free: 400 627 9288

Email: service@genomeditech.com

### **SDS-PAGE**




On SDS-PAGE under reducing (R) condition. The gel was stained overnight with Coomassie Blue. The purity of the protein is greater than 95%.

## **Bioactivity-ELISA**

## **Bioactivity-ELISA**

0.1 μg Human PD1 Protein; His Tag of per well



Human PD1 Protein; His Tag (Catalog # GM-87593RP) was immobilized at 1  $\mu$ g/ml (100  $\mu$ L/well). Increasing concentrations of Anti-PD1 hIgG4 Antibody(Pembrolizumab) (Catalog # GM-52674AB) were added.