Genomeditech (Shanghai) Co.,Ltd.
Order: 021-68455258/50432826/50432825

Toll-free: 400 627 9288

Email: service@genomeditech.com

# Cynomolgus TSLP Protein; His Tag

#### **Product Information**

Product Name Cynomolgus TSLP Protein; His Tag

**Storage temp.** Store at  $\leq$  -70°C, stable for 6 months after receipt.

Recommend to aliquot the protein into smaller quantities for

optimal storage. Please minimize freeze-thaw cycles.

Catalog# / Size GM-87657RP-100 / 100 μg

GM-87657RP-1000 / 1 mg

#### **Protein Information**

Alternative Names Thymic stromal lymphopoietin

Source Cynomolgus TSLP Protein; His Tag (GM-87657RP) is expressed from human

293 cells (HEK-293). It contains AA Tyr 29 - Gln 159 (Accession #

A0A7N9CAT7).

This protein carries a His tag at the C-terminus.

Purity > 95% as determined by SDS-PAGE

Endotoxin < 1 EU/μg, determined by LAL gel clotting assay

Predicted Mol Mass 15.1 KDa

**Formulation** Supplied as a 0.2 μm filtered solution of PBS, pH7.4.

Description TSLP protein (Thymic stromal lymphopoietin) is a cytokine protein mainly

produced by epithelial cells, immune cells, and other cells. It plays an important

role in immune regulation and inflammation processes.

TSLP protein has various functions in the immune system. It can promote the maturation and activation of dendritic cells, thereby triggering immune responses

of T cells and B cells. Additionally, TSLP protein can also affect the activity of

macrophages, mast cells, and other immune cells, regulating inflammatory

responses and immune responses.

TSLP protein plays an important role in allergic diseases and autoimmune

diseases. It is considered an important regulatory factor in allergic diseases (such as asthma, eczema) and autoimmune diseases (such as rheumatoid arthritis,

inflammatory bowel disease). Overactivation of TSLP is often closely related to


the occurrence and development of these diseases.

Genomeditech (Shanghai) Co.,Ltd. Order: 021-68455258/50432826/50432825

Toll-free: 400 627 9288

Email: service@genomeditech.com

#### **SDS-PAGE**




On SDS-PAGE under reducing (R) condition. The gel was stained overnight with Coomassie Blue. The purity of the protein is greater than 95%.

### **Bioactivity-ELISA**

## **Bioactivity-ELISA**

0.1 µg Cynomolgus TSLP Protein; His Tag of per well



Cynomolgus TSLP Protein; His Tag (Catalog # GM-87657RP) was immobilized at 1  $\mu$ g/ml (100  $\mu$ L/well). Increasing concentrations of Anti-TSLP hIgG2 Antibody(Tezepelumab) (Catalog # GM-51675AB) were added.