Genomeditech (Shanghai) Co.,Ltd. Order: 021-68455258/50432826/50432825

Toll-free: 400 627 9288

Email: service@genomeditech.com

# Anti-OX40L hlgG1 Reference Antibody (Oxebio)

## **Product Information**

Product Name Anti-OX40L hlgG1 Reference Antibody (Oxebio)

Storage temp. Store at 2-8°C short term (1-2 weeks). Store at ≤ -20°C long term. Avoid

repeated freeze-thaw.

Catalog# / Size GM-87694MAB-1mg / 1 mg

GM-87694MAB-**5mg / 5 mg** GM-87694MAB-**25mg / 25 mg** GM-87694MAB-**50mg / 50 mg** GM-87694MAB-**100mg / 100 mg** 

# **Antibody Information**

Expression System CHO

Aggregation < 5% as determined by SEC-HPLC

Purity > 95% as determined by SDS-PAGE

Endotoxin < 1 EU/mg, determined by LAL gel clotting assay

Sterility 0.2 µm Filtered

Target OX40L
Clone Oxelumab

Alternative Names CD134L, CD252, GP34, OX-40L, OX4OL, TNLG2B, TXGP1

Source/Isotype Human IgG1(K214R/D356E/L358M, Delta K447), Kappa

Application /

Description OX40L, is not expressed on resting antigen presenting cells, but is

following their activation. Expression of OX40 is dependent on full activation of the T cell; without CD28, expression of OX40 is delayed and of fourfold lower levels. OX40L binds to OX40 receptors on T-cells, preventing them from dying and subsequently increasing cytokine production. OX40 has a critical role in the maintenance of an immune response beyond the first few days and onwards to a memory response due to its ability to enhance survival. OX40 also plays a crucial role in both

Th1 and Th2 mediated reactions in vivo.

**Formulation** phosphate-buffered solution, pH 7.4.

Version: 3.1.3 Revision Date: 15/09/2023


Genomeditech (Shanghai) Co.,Ltd. Order: 021-68455258/50432826/50432825

Toll-free: 400 627 9288

Email: service@genomeditech.com


# **Data Examples**

#### **SDS-PAGE**



On SDS-PAGE under reducing (R)/non-reducing(N-R) condition. The gel was stained overnight with Coomassie Blue. The purity of the protein is greater than 95%.

### **SEC-HPLC**



The purity of this product is more than 95% verified by SEC-HPLC