Genomeditech (Shanghai) Co.,Ltd. Order: 021-68455258/50432826/50432825

Toll-free: 400 627 9288

Email: service@genomeditech.com

Anti-MSLN hlgG1 Reference Antibody (Anetbio)

Product Information

Product Name Anti-MSLN hlgG1 Reference Antibody (Anetbio)

Storage temp. Store at 2-8°C short term (1-2 weeks). Store at ≤ -20°C long term. Avoid

repeated freeze-thaw.

Catalog# / Size GM-87711MAB-1mg / 1 mg

GM-87711MAB-5mg / 5 mg GM-87711MAB-25mg / 25 mg GM-87711MAB-50mg / 50 mg GM-87711MAB-100mg / 100 mg

Antibody Information

Expression System CHO

Aggregation < 5% as determined by SEC-HPLC

Purity > 95% as determined by SDS-PAGE

Endotoxin < 1 EU/mg, determined by LAL gel clotting assay

Sterility 0.2 µm Filtered

Target MSLN

Clone Anetumab

Alternative Names MPF, SMRP

Source/Isotype Human IgG1(K214R), Kappa

Application /

Description MSLN is over expressed in several human tumors, including

mesothelioma, ovarian cancer, pancreatic adenocarcinoma, lung adenocarcinoma, and cholangiocarcinoma. MSLN is a tumor differentiation antigen that is normally present on the mesothelial cells lining the pleura, peritoneum and pericardium. Since mesothelin is overexpressed in several cancers and is immunogenic, the protein could be exploited as tumor marker or as the antigenic target of a therapeutic

cancer vaccine.

Formulation phosphate-buffered solution, pH 7.4.


Genomeditech (Shanghai) Co.,Ltd. Order: 021-68455258/50432826/50432825

Toll-free: 400 627 9288

Email: service@genomeditech.com


Data Examples

SDS-PAGE



On SDS-PAGE under reducing (R)/non-reducing(N-R) condition. The gel was stained overnight with Coomassie Blue. The purity of the protein is greater than 95%.

SEC-HPLC



The purity of this product is more than 95% verified by SEC-HPLC