Genomeditech (Shanghai) Co.,Ltd. Order: 021-68455258/50432826/50432825

Toll-free: 400 627 9288

Email: service@genomeditech.com

# Anti-STEAP2 hlgG1 Reference Antibody (AZD0754)

### **Product Information**

Product Name Anti-STEAP2 hlgG1 Reference Antibody (AZD0754)

**Storage temp.** Store at 2-8°C short term (1-2 weeks). Store at ≤ -20°C long term. Avoid

repeated freeze-thaw.

Catalog# / Size GM-87756MAB-1mg / 1 mg

GM-87756MAB-5mg / 5 mg GM-87756MAB-25mg / 25 mg GM-87756MAB-50mg / 50 mg GM-87756MAB-100mg / 100 mg

# **Antibody Information**

Expression System CHO

Aggregation < 5% as determined by SEC-HPLC

Purity > 95% as determined by SDS-PAGE

Endotoxin < 1 EU/mg, determined by LAL gel clotting assay

Sterility 0.2 µm Filtered

Target STEAP2
Clone AZD0754

Alternative Names IPCA1, PCANAP1, PUMPCn, STAMP1, STMP

Source/Isotype Human IgG1 D356E/L358M, Kappa

Application Flow cytometry

**Description**This gene is a member of the STEAP family and encodes a multi-

pass membrane protein that localizes to the Golgi complex, the plasma membrane, and the vesicular tubular structures in the cytosol. A highly similar protein in mouse has both ferrireductase and cupric reductase activity, and stimulates the cellular uptake of both iron and copper in vitro. Increased transcriptional expression of the human gene is associated with prostate cancer progression. Alternate transcriptional splice variants, encoding different

isoforms, have been characterized.

**Formulation** phosphate-buffered solution, pH 7.4.


Genomeditech (Shanghai) Co.,Ltd. Order: 021-68455258/50432826/50432825

Toll-free: 400 627 9288

Email: service@genomeditech.com


# **Data Examples**

# SDS-PAGE



On SDS-PAGE under reducing (R)/non-reducing(N-R) condition. The gel was stained overnight with Coomassie Blue. The purity of the protein is greater than 95%.

# **SEC-HPLC**



The purity of this product is more than 95% verified by SEC-HPLC




Genomeditech (Shanghai) Co.,Ltd. Order: 021-68455258/50432826/50432825

Toll-free: 400 627 9288

Email: service@genomeditech.com

Flow cytometry

H\_STEAP2(ECD) CHO-K1 Cell Line (Catalog # GM-C37055) was stained with Anti-STEAP2 hlgG1 Reference Antibody (AZD0754) (Catalog # GM-87756MAB) or isotype control antibody, followed by anti-Human IgG APC-conjugated Secondary Antibody.



SampleID	Geometric Mean : FL11-H
CHO-K1 anti-STEAP2+APC-2nd Ab	2340
CHO-K1 H_STEAP2(ECD) H_IgG+APC-2nd Ab	2546
CHO-K1 H_STEAP2(ECD) anti-STEAP2 +APC-2nd Ab	85111

Fig. FACS