Genomeditech (Shanghai) Co.,Ltd. Order: 021-68455258/50432826/50432825

Toll-free: 400 627 9288

Email: service@genomeditech.com

Anti-PD-L1 hlgG1 Reference Antibody (Avebio)

Product Information

Product Name Anti-PD-L1 hlgG1 Reference Antibody (Avebio)

Storage temp. Store at $2-8^{\circ}$ C short term (1-2 weeks). Store at $\leq -20^{\circ}$ C long term. Avoid

repeated freeze-thaw.

Catalog# / Size GM-87782MAB-1mg / 1 mg

GM-87782MAB-5mg / 5 mg GM-87782MAB-25mg / 25 mg GM-87782MAB-50mg / 50 mg GM-87782MAB-100mg / 100 mg

Antibody Information

Expression System CHO

Aggregation < 5% as determined by SEC-HPLC

Purity > 95% as determined by SDS-PAGE

Endotoxin < 1 EU/mg, determined by LAL gel clotting assay

Sterility 0.2 µm Filtered

Target PD-L1
Clone Avelumab

Alternative Names B7-H, B7H1, PDL1, PDCD1L1, PDCD1LG1, CD274, hPD-L1

Source/Isotype Human IgG1 (KDEL), Kappa

Application Bioactivity-ELISA

Description The programmed cell death 1 protein (PD-1, PDCD1, CD279) is a member

of the CD28 family of immunoreceptors that regulate T cell activation and immune responses. The PD-1 protein contains an extracellular Ig V domain, a transmembrane domain, and a cytoplasmic tail that includes an immunoreceptor tyrosine-based inhibitory motif (ITIM) and an immunoreceptor tyrosine-based switch motif (ITSM). PD-1 is activated by the cell surface ligands PD-L1 and PD-L2. Upon activation, PD-1 ITIM and ITSM phosphorylation leads to the recruitment of the protein tyrosine phosphatases SHP-1 and SHP-2, which suppress TCR signaling. In addition to activated T-cells, PD-1 is expressed in activated B-cells and monocytes, although its function in these cell types has not been fully characterized. The PD-1 pathway plays an important role in immune tolerance; however, research studies show that cancer cells often adopt this pathway to escape immune surveillance. Consequently, blockade of

Version:3.2



Genomeditech (Shanghai) Co.,Ltd. Order: 021-68455258/50432826/50432825

Toll-free: 400 627 9288

Email: service@genomeditech.com

PD-1 and its ligands is proving to be a sound strategy for neoplastic


intervention.

Formulation

phosphate-buffered solution, pH 7.4.


Data Examples

SDS-PAGE



On SDS-PAGE under reducing (R)/non-reducing(N-R) condition. The gel was stained overnight with Coomassie Blue. The purity of the protein is greater than 95%.

SEC-HPLC



The purity of this product is more than 95% verified by SEC-HPLC

Genomeditech (Shanghai) Co.,Ltd. Order: 021-68455258/50432826/50432825

Toll-free: 400 627 9288


Email: service@genomeditech.com

Bioactivity-ELISA

Human PDL1 Protein; His Tag (Catalog # GM-85609RP) was immobilized at 1 μ g/ml (100 μ L/well). Increasing concentrations of Anti-PD-L1 hlgG1 Reference Antibody (Avebio) (Catalog # GM-87782MAB) were added.

Bioactivity-ELISA

0.1 µg Human PDL1 Protein; His Tag of per well



Anti-PD-L1 hIgG1 Reference Antibody(Avebio).(ng/mL)

Fig. Assay

Bioactivity-ELISA

Cynomolgus PDL1 Protein; His Tag (Catalog # GM-85611RP) was immobilized at 1 μ g/ml (100 μ L/well) on Anti-His mlgG2a Antibody (Catalog # GM-59493AB) (0.6 μ g/well) precoated. Increasing concentrations of Anti-PD-L1 hlgG1 Reference Antibody (Avebio) (Catalog # GM-87782MAB) were added.

Bioactivity-ELISA

0.6 μg Anti-His mIgG2a Antibody+0.1 μg Cynomolgus PDL1 Protein; His Tag of per well




Fig. Assay