Genomeditech (Shanghai) Co.,Ltd. Order: 021-68455258/50432826/50432825

Toll-free: 400 627 9288

Email: service@genomeditech.com

Anti-TL1A hlgG1 Reference Antibody (Duvbio)

Product Information

Product Name Anti-TL1A hlgG1 Reference Antibody (Duvbio)

Storage temp. Store at $2-8^{\circ}$ C short term (1-2 weeks). Store at $\leq -20^{\circ}$ C long term. Avoid

repeated freeze-thaw.

Catalog# / Size GM-87806MAB-1mg / 1 mg

GM-87806MAB-5mg / 5 mg GM-87806MAB-25mg / 25 mg GM-87806MAB-50mg / 50 mg GM-87806MAB-100mg / 100 mg

Antibody Information

Expression System CHO

Aggregation < 5% as determined by SEC-HPLC

Purity > 95% as determined by SDS-PAGE

Endotoxin < 1 EU/mg, determined by LAL gel clotting assay

Sterility 0.2 µm Filtered

Target TL1A
Clone Duvakitug

Alternative Names TL1; TL1A; VEGI; TNLG1B; VEGI192A

Source/Isotype Human IgG1 (KDEL), Lamda2

Application /

Description The protein encoded by TNFSF15 is a cytokine that belongs to the tumor

necrosis factor (TNF) ligand family. This protein is abundantly expressed in endothelial cells but is not expressed in either B or T cells. The expression of this protein is inducible by TNF-alpha and IL-1 alpha. This cytokine is a ligand for receptor TNFRSF25 (death receptor 3) and TNFRSF6B (decoy receptor 3). It can activate both the NF-kB and MAPK signalling pathways and acts as an autocrine factor to induce apoptosis in endothelial cells. This cytokine is also found to inhibit endothelial cell proliferation, and thus may function as an angiogenesis inhibitor. An additional isoform encoded by an alternatively spliced transcript variant has been reported but the sequence of this transcript has not been

determined.

Formulation phosphate-buffered solution, pH 7.4.

Version:3.1


Genomeditech (Shanghai) Co.,Ltd. Order: 021-68455258/50432826/50432825

Toll-free: 400 627 9288

Email: service@genomeditech.com


Data Examples

SDS-PAGE



On SDS-PAGE under reducing (R)/non-reducing(N-R) condition. The gel was stained overnight with Coomassie Blue. The purity of the protein is greater than 95%.

SEC-HPLC



The purity of this product is more than 95% verified by SEC-HPLC