

Anti-CD40 hlgG1 Reference Antibody (Tenebio)

Product Information


Product Name	Anti-CD40 hlgG1 Reference Antibody (Tenebio)
Storage temp.	Store at 2-8°C short term (1-2 weeks).Store at $\leq -20^{\circ}\text{C}$ long term. Avoid repeated freeze-thaw.
Catalog# / Size	GM-87983MAB-1mg / 1 mg GM-87983MAB-5mg / 5 mg GM-87983MAB-25mg / 5 mg * 5 vials GM-87983MAB-50mg / 50 mg GM-87983MAB-100mg / 50 mg * 2 vials

Antibody Information

Expression System	CHO
Aggregation	< 5% as determined by SEC-HPLC
Purity	> 95% as determined by SDS-PAGE
Endotoxin	< 1 EU/mg, determined by LAL gel clotting assay
Sterility	0.2 μm Filtered
Target	CD40
Clone	Teneliximab
Alternative Names	Bp50, CDW40, TNFRSF5, p50
Source/Isotype	Monoclonal Human IgG1(KDEL), Kappa
Application	/
Description	CD40 (also known as Bp50, CDW40, TNFRSF5, and p50) is a membrane protein receptor encoded by the CD40 gene, located on human chromosome 20q13. Teneliximab is a monoclonal antibody targeting CD40, primarily functioning to regulate the CD40/CD40L signaling pathway. By blocking or modulating CD40 activity, Teneliximab plays a significant role in cancer immunotherapy and the treatment of immune-related diseases.
Formulation	phosphate-buffered solution, pH 7.2-7.4.


Data Examples

SDS-PAGE



On SDS-PAGE under reducing (R)/non-reducing(N-R) condition. The gel was stained overnight with Coomassie Blue. The purity of the protein is greater than 95%.

SEC-HPLC



The purity of this product is more than 95% verified by SEC-HPLC