

Human CDCP1 Protein; hFc Tag


Product Information

Product Name	Human CDCP1 Protein; hFc Tag
Storage temp	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.
Catalog# / Size	GM-88172RP-100 / 100 μg GM-88172RP-1000 / 1 mg

Protein Information


Alternative Names	CDCP1, CD318, SIMA135, TRASK, UNQ2486, PRO5773
Source	Human CDCP1 Protein; hFc Tag (GM-88172RP) is expressed from human 293 cells (HEK-293). It contains AA Phe 30 - Thr 667 (Accession # Q9H5V8-1). This protein carries a hFc tag at the C-terminus.
Purity	> 95% as determined by SDS-PAGE
Endotoxin	< 1 EU/ μg , determined by LAL gel clotting assay
Predicted Mol Mass	97.8 KDa
Formulation	Supplied as a 0.2 μm filtered solution of PBS, pH7.2-7.4.
Description	CDCP1 (CD318) is a transmembrane glycoprotein broadly expressed on epithelial cells in skin, adrenal glands, and intestines, plus some blood cells and tumors. It features a long extracellular domain, a single transmembrane segment, and a short cytoplasmic tail. CDCP1 participates in adhesion, migration, signaling, and remodeling, activating pathways like PI3K/Akt and Src to influence survival and movement. Clinically, CDCP1 links to cancer progression, especially invasion and metastasis. Its level and activation may serve as a prognostic marker or therapeutic target, possibly by modulating cell-cell and cell-matrix interactions and the tumor microenvironment. Research continues on its upstream regulators, downstream signaling, and therapies such as antibodies, inhibitors, or antibody-drug conjugates to block its functions.

SDS-PAGE



On SDS-PAGE under non-reducing (NR) condition and reducing (R) condition. The gel was stained overnight with Coomassie Blue. The purity of the protein is greater than 95%.

Bioactivity-ELISA



Human CDCP1 Protein; hFc Tag (Catalog # GM-88172RP) was immobilized at 2 µg/ml (100 µL/well). Increasing concentrations of Anti-CDCP1 mIgG2a Antibody (38E11) (Catalog # GM-88178AB) were added.